

Product Name	Description	Size	Price	Intended Use
VS208070(A)	Western blot stripping buffer (1x)	x 2	\$280.00	Effectively remove all 1st and 2nd antibodies on membranes (PVDF/NC)
VS208070	Western blot stripping buffer (1x)	500 mL	\$120.00	Effectively remove all 1st and 2nd antibodies on membranes (PVDF/NC)
VR208090	Protein Signal Recover Buffer	125 mL 10x	\$150.00	Remove white bands/yellow bands/empty, uneven bands on x-films
VR208090(A)	Protein Signal Recover Buffer	500 mL (1x)	\$120.00	Remove white bands/yellow bands/empty, uneven bands on x-films
VF208010	X-ray Film background remover reagents (A:B)	250 mL	\$250.00	Remove backgrounds on x-ray films; overcome the film over-exposed & saturation
C218060	Western blot background reducing buffer	500 mL	\$150.00	Remove high backgrounds on membranes. 5-10 minutes washing before ECL development
VB208068	Western blot blocking buffer (New)	500 mL	\$150.00	Block no-specific backgrounds, enhance antibody bindings, and reduce high backgrounds
VA212060	Affinity Enhancer Buffer (A:B)	500 mL	\$150.00	Effectively enhance 1st antibody binding and reduce no-specific bindings